

Recombinant Human GRP Protein, N-His-KSI

Catalog No: #AP98135

Package Size: #AP98135-1 50ug #AP98135-2 100ug #AP98135-3 1mg

Description

Product Name	Recombinant Human GRP Protein, N-His-KSI
Host Species	E. coli
Purification	>90% as determined by SDS-PAGE.
Applications	ELISA, Immunogen, SDS-PAGE, WB, Bioactivity testing in progress
Other Names	GRP-10, Gastrin-releasing peptide, GRP
Accession No.	P07492
Target Species	Homo sapiens (Human)
Calculated MW	29.47 kDa
Target Length	Val24-Gln148
Tag Info	His-tag
Formulation	Lyophilized from a solution in PBS pH 7.4, 0.02% NLS, 1mM EDTA, 4% Trehalose, 1% Mannitol.
Storage	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for one week. Store at -20 to -80°C for twelve months from the date of receipt.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.