

Recombinant Human MPZ Protein, N-His

Catalog No: #AP98137

Package Size: #AP98137-1 50ug #AP98137-2 100ug #AP98137-3 1mg

Description

Product Name	Recombinant Human MPZ Protein, N-His
Host Species	E. coli
Purification	>90% as determined by SDS-PAGE.
Applications	ELISA, Immunogen, SDS-PAGE, WB, Bioactivity testing in progress
Other Names	MPZ, Myelin protein zero, Myelin protein P0, Myelin peripheral protein, MPP
Accession No.	P25189
Target Species	Homo sapiens (Human)
Calculated MW	15.93 kDa
Target Length	Ile30-Glu148
Tag Info	His-tag
Formulation	Lyophilized from a solution in PBS pH 7.4, 0.02% NLS, 1mM EDTA, 4% Trehalose, 1% Mannitol.
Storage	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for one week. Store at -20 to -80°C for twelve months from the date of receipt.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.