

Recombinant Human S100A4 Protein, N-GST & C-His

Catalog No: #AP98138

Package Size: #AP98138-1 50ug #AP98138-2 100ug #AP98138-3 1mg

Description

Product Name	Recombinant Human S100A4 Protein, N-GST & C-His
Host Species	E. coli
Purification	>90% as determined by SDS-PAGE.
Applications	ELISA, Immunogen, SDS-PAGE, WB, Bioactivity testing in progress
Other Names	S100A4, Protein S100-A4, S100 calcium-binding protein A4, Calvasculin, Protein Mts1, MTS1, CAPL, Placental calcium-binding protein, Metastasin
Accession No.	P26447
Target Species	Homo sapiens (Human)
Calculated MW	39.78 kDa
Target Length	Ala2-Lys101
Formulation	Lyophilized from a solution in PBS pH 7.4, 0.02% NLS, 1mM EDTA, 4% Trehalose, 1% Mannitol.
Storage	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for one week. Store at -20 to -80°C for twelve months from the date of receipt.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.