

Recombinant Human SELENOW Protein, N-GST & C-His

Catalog No: #AP98142

Package Size: #AP98142-1 50ug #AP98142-2 100ug #AP98142-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Human SELENOW Protein, N-GST & C-His
Host Species	E. coli
Purification	>90% as determined by SDS-PAGE.
Applications	ELISA, Immunogen, SDS-PAGE, WB, Bioactivity testing in progress
Other Names	SELW, Selenoprotein W, SEPW1, SELENOW, SelW
Accession No.	P63302
Target Species	Homo sapiens (Human)
Calculated MW	37.58 kDa
Target Length	Met1-Gly87
Formulation	Lyophilized from a solution in PBS pH 7.4, 0.02% NLS, 1mM EDTA, 4% Trehalose, 1% Mannitol.
Storage	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for one week. Store at -20 to -80°C for twelve months from the date of receipt.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.