

Recombinant Human CCL19/MIP-3-beta Protein, N-GST & C-His

Catalog No: #AP98145

Package Size: #AP98145-1 50ug #AP98145-2 100ug #AP98145-3 1mg

Description

Product Name	Recombinant Human CCL19/MIP-3-beta Protein, N-GST & C-His
Host Species	E. coli
Purification	>90% as determined by SDS-PAGE.
Applications	ELISA, Immunogen, SDS-PAGE, WB, Bioactivity testing in progress
Other Names	C-C motif chemokine 19, MIP-3-beta, Macrophage inflammatory protein 3 beta, Small-inducible cytokine A19, EBI1 ligand chemokine, SCYA19, CK beta-11, Beta-chemokine exodus-3, Epstein-Barr virus-induced molecule 1 ligand chemokine, CCL19, MIP3B, ELC
Accession No.	Q99731
Target Species	Homo sapiens (Human)
Calculated MW	36.98 kDa
Target Length	Gly22-Ser98
Formulation	Lyophilized from a solution in PBS pH 7.4, 0.02% NLS, 1mM EDTA, 4% Trehalose, 1% Mannitol.
Storage	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for one week. Store at -20 to -80°C for twelve months from the date of receipt.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.