

Recombinant HSV-1/HHV-1 gB Protein, C-Strep

Catalog No: #AP98148

Package Size: #AP98148-1 50ug #AP98148-2 100ug #AP98148-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant HSV-1/HHV-1 gB Protein, C-Strep
Host Species	Mammalian cells
Purification	>90% as determined by SDS-PAGE.
Applications	ELISA, Immunogen, SDS-PAGE, WB, Bioactivity testing in progress
Other Names	gB, Envelope glycoprotein B
Accession No.	A0A181ZGK4
Target Species	Human herpesvirus 1 (strain 17) (HHV-1) (Human herpes simplex virus 1)
Calculated MW	81.65 kDa
Target Length	Ala30-Ala730
Tag Info	Strep-tag
Formulation	Lyophilized from a solution in PBS pH 7.4, 1mM EDTA, 4% Trehalose, 1% Mannitol.
Storage	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for one week. Store at -20 to -80°C for twelve months from the date of receipt.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.