

Recombinant HSV-1/HHV-1 gE Protein, C-His

Catalog No: #AP98151

Package Size: #AP98151-1 50ug #AP98151-2 100ug #AP98151-3 1mg

Description

Product Name	Recombinant HSV-1/HHV-1 gE Protein, C-His
Host Species	Mammalian cells
Purification	>90% as determined by SDS-PAGE.
Applications	ELISA, Immunogen, SDS-PAGE, WB, Bioactivity testing in progress
Other Names	gE, Envelope glycoprotein E
Accession No.	P04488
Target Species	Human herpesvirus 1 (strain 17) (HHV-1) (Human herpes simplex virus 1)
Calculated MW	45.76 kDa
Target Length	Gly21-Arg419
Tag Info	His-tag
Formulation	Lyophilized from a solution in PBS pH 7.4, 1mM EDTA, 4% Trehalose, 1% Mannitol.
Storage	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for one week. Store at -20 to -80°C for twelve months from the date of receipt.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.