

Recombinant HSV-1/HHV-1 gG/US4 Protein, N-GST

Catalog No: #AP98155

Package Size: #AP98155-1 50ug #AP98155-2 100ug #AP98155-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant HSV-1/HHV-1 gG/US4 Protein, N-GST
Host Species	E. coli
Purification	>90% as determined by SDS-PAGE.
Applications	ELISA, Immunogen, SDS-PAGE, WB, Bioactivity testing in progress
Other Names	gG, US4
Accession No.	P06484
Target Species	Human herpesvirus 1 (strain 17) (HHV-1) (Human herpes simplex virus 1)
Calculated MW	35.98 kDa
Target Length	Glu84-His175
Tag Info	GST-tag
Formulation	Lyophilized from a solution in PBS pH 7.4, 0.02% NLS, 1mM EDTA, 4% Trehalose, 1% Mannitol.
Storage	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for one week. Store at -20 to -80°C for twelve months from the date of receipt.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.