

Recombinant HSV-2/HHV-2 gC/UL44 Protein, C-His

Catalog No: #AP98159

Package Size: #AP98159-1 50ug #AP98159-2 100ug #AP98159-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant HSV-2/HHV-2 gC/UL44 Protein, C-His
Host Species	Mammalian cells
Purification	>90% as determined by SDS-PAGE.
Applications	ELISA, Immunogen, SDS-PAGE, WB, Bioactivity testing in progress
Other Names	gC, UL44, Envelope glycoprotein C
Accession No.	A0A0E3Y7C2
Target Species	Human herpesvirus 2 (HHV-2) (Human herpes simplex virus 2)
Calculated MW	48.43 kDa
Target Length	Ala27-Gly447
Tag Info	His-tag
Formulation	Lyophilized from a solution in PBS pH 7.4, 1mM EDTA, 4% Trehalose, 1% Mannitol.
Storage	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for one week. Store at -20 to -80°C for twelve months from the date of receipt.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.