

## Recombinant JUNV GP1 Protein, C-Fc

Catalog No: #AP98162

Package Size: #AP98162-1 50ug #AP98162-2 100ug #AP98162-3 1mg

## Description

|                |                                                                                                                                                                       |
|----------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Product Name   | Recombinant JUNV GP1 Protein, C-Fc                                                                                                                                    |
| Host Species   | Mammalian cells                                                                                                                                                       |
| Purification   | >90% as determined by SDS-PAGE.                                                                                                                                       |
| Applications   | ELISA, Immunogen, SDS-PAGE, WB, Bioactivity testing in progress                                                                                                       |
| Other Names    | Pre-glycoprotein polyprotein GP complex, Glycoprotein G1, GP1, GPC, GP-C                                                                                              |
| Accession No.  | P26313                                                                                                                                                                |
| Target Species | Junin mammarenavirus (JUNV) (Junn mammarenavirus)                                                                                                                     |
| Calculated MW  | 50.49 kDa                                                                                                                                                             |
| Target Length  | Glu59-Leu250                                                                                                                                                          |
| Tag Info       | Fc-tag                                                                                                                                                                |
| Formulation    | Lyophilized from a solution in PBS pH 7.4, 1mM EDTA, 4% Trehalose, 1% Mannitol.                                                                                       |
| Storage        | Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for one week. Store at -20 to -80°C for twelve months from the date of receipt. |

Note: This product is for in vitro research use only and is not intended for use in humans or animals.