

Recombinant JUNV GP2 Protein, C-Fc

Catalog No: #AP98163

Package Size: #AP98163-1 50ug #AP98163-2 100ug #AP98163-3 1mg

Description

Product Name	Recombinant JUNV GP2 Protein, C-Fc
Host Species	Mammalian cells
Purification	>90% as determined by SDS-PAGE.
Applications	ELISA, Immunogen, SDS-PAGE, WB, Bioactivity testing in progress
Other Names	Pre-glycoprotein polyprotein GP complex, Glycoprotein G2, GP2, GPC, GP-C
Accession No.	P26313
Target Species	Junin mammarenavirus (JUNV) (Junn mammarenavirus)
Calculated MW	48.36 kDa
Target Length	Ala252-Asp424
Tag Info	Fc-tag
Formulation	Lyophilized from a solution in PBS pH 7.4, 1mM EDTA, 4% Trehalose, 1% Mannitol.
Storage	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for one week. Store at -20 to -80°C for twelve months from the date of receipt.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.