

Recombinant Influenza A virus (H3N2) NR (N)-3M2e Protein, N-His

Catalog No: #AP98176

Package Size: #AP98176-1 50ug #AP98176-2 100ug #AP98176-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Influenza A virus (H3N2) NR (N)-3M2e Protein, N-His
Host Species	E. coli
Purification	>90% as determined by SDS-PAGE.
Applications	ELISA, Immunogen, SDS-PAGE, WB, Bioactivity testing in progress
Other Names	Matrix protein 2, Proton channel protein M2, M, Nucleoprotein, Protein N, Nucleocapsid protein, N
Accession No.	P03418 & BFM37795.1
Target Species	Influenza A virus H3N2
Calculated MW	54.72 kDa
Target Length	[Met1-Leu391] & [Met1-Pro25]3
Tag Info	His-tag
Formulation	Lyophilized from a solution in PBS pH 7.4, 0.02% NLS, 1mM EDTA, 4% Trehalose, 1% Mannitol.
Storage	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for one week. Store at -20 to -80°C for twelve months from the date of receipt.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.