

Recombinant Mouse CD132/IL2RG Protein, N-His

Catalog No: #AP98186

Package Size: #AP98186-1 50ug #AP98186-2 100ug #AP98186-3 1mg

Description

Product Name	Recombinant Mouse CD132/IL2RG Protein, N-His
Host Species	E. coli
Purification	>90% as determined by SDS-PAGE.
Applications	ELISA, Immunogen, SDS-PAGE, WB, Bioactivity testing in progress
Other Names	Cytokine receptor common subunit gamma, Interleukin-2 receptor subunit gamma, IL-2 receptor subunit gamma, IL-2R subunit gamma, IL-2RG, gammaC, p64, CD132, Il2rg
Accession No.	P34902
Target Species	Mus musculus (Mouse)
Calculated MW	30.52 kDa
Target Length	Trp23-Ala263
Tag Info	His-tag
Formulation	Lyophilized from a solution in PBS pH 7.4, 0.02% NLS, 1mM EDTA, 4% Trehalose, 1% Mannitol.
Storage	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for one week. Store at -20 to -80°C for twelve months from the date of receipt.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.