

Recombinant Treponema pallidum FlaB3 Protein, N-His-SUMO

Catalog No: #AP98188

Package Size: #AP98188-1 50ug #AP98188-2 100ug #AP98188-3 1mg

Description

Product Name	Recombinant Treponema pallidum FlaB3 Protein, N-His-SUMO
Host Species	E. coli
Purification	>90% as determined by SDS-PAGE.
Applications	ELISA, Immunogen, SDS-PAGE, WB, Bioactivity testing in progress
Other Names	Flagellin FlaB3, Class B, Flagellar filament 31 kDa core protein, flaB3, TP_0870
Accession No.	P21992
Target Species	Treponema pallidum (strain Nichols)
Calculated MW	43.41 kDa
Target Length	Met1-Gln285
Tag Info	His-tag
Formulation	Lyophilized from a solution in PBS pH 7.4, 0.02% NLS, 1mM EDTA, 4% Trehalose, 1% Mannitol.
Storage	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for one week. Store at -20 to -80°C for twelve months from the date of receipt.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.