

Recombinant Human CCKBR/CCK2-R Protein, N-His

Catalog No: #AP98197

Package Size: #AP98197-1 50ug #AP98197-2 100ug #AP98197-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Human CCKBR/CCK2-R Protein, N-His
Host Species	E. coli
Purification	>90% as determined by SDS-PAGE.
Applications	ELISA, Immunogen, SDS-PAGE, WB, Bioactivity testing in progress
Other Names	Cholecystokinin-2 receptor, CCK2-R, CCK-BR, CCK-B receptor, Gastrin/cholecystokinin type B receptor, CCKBR, CCKRB
Accession No.	P32239
Target Species	Homo sapiens (Human)
Calculated MW	12.00 kDa
Target Length	Arg243-Arg333
Tag Info	His-tag
Formulation	Lyophilized from a solution in PBS pH 7.4, 0.02% NLS, 1mM EDTA, 4% Trehalose, 1% Mannitol.
Storage	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for one week. Store at -20 to -80°C for twelve months from the date of receipt.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.