

Recombinant EBV/HHV-4 gB/BALF4 Protein, N-His

Catalog No: #AP98200

Package Size: #AP98200-1 50ug #AP98200-2 100ug #AP98200-3 1mg

Description

Product Name	Recombinant EBV/HHV-4 gB/BALF4 Protein, N-His
Host Species	E. coli
Purification	>90% as determined by SDS-PAGE.
Applications	ELISA, Immunogen, SDS-PAGE, WB, Bioactivity testing in progress
Other Names	Envelope glycoprotein B, gB, BALF4, gB, EBVaGC1_061
Accession No.	R4R670
Target Species	Epstein-Barr virus (strain GD1) (HHV-4) (Human herpesvirus 4)
Calculated MW	26.24 kDa
Target Length	Ile90-Cys295 (Glu128Asp, Thr241Lys)
Tag Info	His-tag
Formulation	Lyophilized from a solution in PBS pH 7.4, 0.02% NLS, 1mM EDTA, 4% Trehalose, 1% Mannitol.
Storage	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for one week. Store at -20 to -80°C for twelve months from the date of receipt.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.