

Recombinant Human BTNL2/BTL-II Protein, C-Strep

Catalog No: #AP98204

Package Size: #AP98204-1 50ug #AP98204-2 100ug #AP98204-3 1mg

Description

Product Name	Recombinant Human BTNL2/BTL-II Protein, C-Strep
Host Species	Insect Cells
Purification	>90% as determined by SDS-PAGE.
Applications	ELISA, Immunogen, SDS-PAGE, WB, Bioactivity testing in progress
Other Names	BTNL2, Butyrophilin-like protein 2, BTL-II
Accession No.	Q9UIR0
Target Species	Homo sapiens (Human)
Calculated MW	50.51 kDa
Target Length	Phe29-Ser455
Tag Info	Strep-tag
Formulation	Lyophilized from a solution in PBS pH 7.4, 1mM EDTA, 4% Trehalose, 1% Mannitol.
Storage	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for one week. Store at -20 to -80°C for twelve months from the date of receipt.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.