

Recombinant Human CD35/CR1 Protein, C-Fc

Catalog No: #AP98205

Package Size: #AP98205-1 50ug #AP98205-2 100ug #AP98205-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Human CD35/CR1 Protein, C-Fc
Host Species	Mammalian cells
Purification	>90% as determined by SDS-PAGE.
Applications	ELISA, Immunogen, SDS-PAGE, WB, Bioactivity testing in progress
Other Names	C3BR, Complement receptor type 1, CR1, CD35, C3b/C4b receptor
Accession No.	P17927
Target Species	Homo sapiens (Human)
Calculated MW	49.93 kDa
Target Length	Gly941-Asn1136
Tag Info	Fc-tag
Formulation	Lyophilized from a solution in PBS pH 7.4, 1mM EDTA, 4% Trehalose, 1% Mannitol.
Storage	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for one week. Store at -20 to -80°C for twelve months from the date of receipt.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.