

Recombinant Tupaia chinensis CD366/HAVCR2/TIM-3 Protein, C-Strep

Catalog No: #AP98207

Package Size: #AP98207-1 50ug #AP98207-2 100ug #AP98207-3 1mg

Description

Product Name	Recombinant Tupaia chinensis CD366/HAVCR2/TIM-3 Protein, C-Strep
Host Species	Mammalian cells
Purification	>90% as determined by SDS-PAGE.
Applications	ELISA, Immunogen, SDS-PAGE, WB, Bioactivity testing in progress
Other Names	T-cell immunoglobulin mucin receptor 3, HAVcr-2, TIMD-3, T-cell membrane protein 3, HAVCR2, Hepatitis A virus cellular receptor 2, CD366, TIMD3, TIM-3, TIM3, T-cell immunoglobulin and mucin domain-containing protein 3,
Accession No.	XP_014438264.1
Target Species	Tupaia chinensis (Chinese tree shrew)
Calculated MW	21.99 kDa
Target Length	Leu23-Arg194
Tag Info	Strep-tag
Formulation	Lyophilized from a solution in PBS pH 7.4, 1mM EDTA, 4% Trehalose, 1% Mannitol.
Storage	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for one week. Store at -20 to -80°C for twelve months from the date of receipt.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.