

Recombinant Tupaia chinensis CD137/TNFRSF9/4-1BB Protein, C-Fc

Catalog No: #AP98214

Package Size: #AP98214-1 50ug #AP98214-2 100ug #AP98214-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Tupaia chinensis CD137/TNFRSF9/4-1BB Protein, C-Fc
Host Species	Mammalian cells
Purification	>90% as determined by SDS-PAGE.
Applications	ELISA, Immunogen, SDS-PAGE, WB, Bioactivity testing in progress
Other Names	Tumor necrosis factor receptor superfamily member 9, CD137, 4-1BB ligand receptor, CDw137, TNFRSF9, T-cell antigen ILA, T-cell antigen 4-1BB homolog, ILA,
Accession No.	XP_006145444.1
Target Species	Tupaia chinensis (Chinese tree shrew)
Calculated MW	44.53 kDa
Target Length	Thr29-Thr180
Tag Info	Fc-tag
Formulation	Lyophilized from a solution in PBS pH 7.4, 1mM EDTA, 4% Trehalose, 1% Mannitol.
Storage	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for one week. Store at -20 to -80°C for twelve months from the date of receipt.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.