

Recombinant Tupaia chinensis CD274/PD-L1/B7-H1 Protein, C-Fc

Catalog No: #AP98220

Package Size: #AP98220-1 50ug #AP98220-2 100ug #AP98220-3 1mg

Description

Product Name	Recombinant Tupaia chinensis CD274/PD-L1/B7-H1 Protein, C-Fc
Host Species	Mammalian cells
Purification	>90% as determined by SDS-PAGE.
Applications	ELISA, Immunogen, SDS-PAGE, WB, Bioactivity testing in progress
Other Names	B7-H1, Programmed cell death 1 ligand 1, PDCD1 ligand 1, PDCD1L1, B7 homolog 1, PDCD1LG1, PDL1, hPD-L1, Programmed death ligand 1, B7H1, PD-L1, CD274,
Accession No.	XP_014444807.1
Target Species	Tupaia chinensis (Chinese tree shrew)
Calculated MW	53.08 kDa
Target Length	Phe19-Arg236
Tag Info	Fc-tag
Formulation	Lyophilized from a solution in PBS pH 7.4, 1mM EDTA, 4% Trehalose, 1% Mannitol.
Storage	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for one week. Store at -20 to -80°C for twelve months from the date of receipt.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.