

# Recombinant Tupaia chinensis CD47/MER6 Protein, C-Fc

Catalog No: #AP98221

Package Size: #AP98221-1 50ug #AP98221-2 100ug #AP98221-3 1mg

## Description

|                |                                                                                                                                                                       |
|----------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Product Name   | Recombinant Tupaia chinensis CD47/MER6 Protein, C-Fc                                                                                                                  |
| Host Species   | Mammalian cells                                                                                                                                                       |
| Purification   | >90% as determined by SDS-PAGE.                                                                                                                                       |
| Applications   | ELISA, Immunogen, SDS-PAGE, WB, Bioactivity testing in progress                                                                                                       |
| Other Names    | Leukocyte surface antigen CD47, CD47, Antigenic surface determinant protein OA3, IAP, Integrin-associated protein, Protein MER6, MER6,                                |
| Accession No.  | XP_027625004.1                                                                                                                                                        |
| Target Species | Tupaia chinensis (Chinese tree shrew)                                                                                                                                 |
| Calculated MW  | 41.93 kDa                                                                                                                                                             |
| Target Length  | Gln373-Glu493                                                                                                                                                         |
| Tag Info       | Fc-tag                                                                                                                                                                |
| Formulation    | Lyophilized from a solution in PBS pH 7.4, 1mM EDTA, 4% Trehalose, 1% Mannitol.                                                                                       |
| Storage        | Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for one week. Store at -20 to -80°C for twelve months from the date of receipt. |

Note: This product is for in vitro research use only and is not intended for use in humans or animals.