

Recombinant Tupaia chinensis CD8A Protein, N-His

Catalog No: #AP98223

Package Size: #AP98223-1 50ug #AP98223-2 100ug #AP98223-3 1mg

Description

Product Name	Recombinant Tupaia chinensis CD8A Protein, N-His
Host Species	E. coli
Purification	>90% as determined by SDS-PAGE.
Applications	ELISA, Immunogen, SDS-PAGE, WB, Bioactivity testing in progress
Other Names	T-cell surface glycoprotein CD8 alpha chain, CD8a, T-lymphocyte differentiation antigen T8/Leu-2, MAL, CD8A,
Accession No.	XP_014439628.1
Target Species	Tupaia chinensis (Chinese tree shrew)
Calculated MW	16.25 kDa
Target Length	Trp1-Asp131
Tag Info	His-tag
Formulation	Lyophilized from a solution in PBS pH 7.4, 0.02% NLS, 1mM EDTA, 4% Trehalose, 1% Mannitol.
Storage	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for one week. Store at -20 to -80°C for twelve months from the date of receipt.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.