

Recombinant Tupaia chinensis CD4 Protein, C-His

Catalog No: #AP98225

Package Size: #AP98225-1 50ug #AP98225-2 100ug #AP98225-3 1mg

Description

Product Name	Recombinant Tupaia chinensis CD4 Protein, C-His
Host Species	Mammalian cells
Purification	>90% as determined by SDS-PAGE.
Applications	ELISA, Immunogen, SDS-PAGE, WB, Bioactivity testing in progress
Other Names	T-cell surface antigen T4/Leu-3, T-cell surface glycoprotein CD4, CD4,
Accession No.	XP_014440610.1
Target Species	Tupaia chinensis (Chinese tree shrew)
Calculated MW	43.96 kDa
Target Length	Lys26-Thr394
Tag Info	His-tag
Formulation	Lyophilized from a solution in PBS pH 7.4, 1mM EDTA, 4% Trehalose, 1% Mannitol.
Storage	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for one week. Store at -20 to -80°C for twelve months from the date of receipt.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.