

Recombinant Mouse F2RL1/PAR-2 Protein, N-GST & C-His

Catalog No: #AP98231

Package Size: #AP98231-1 50ug #AP98231-2 100ug #AP98231-3 1mg

Description

Product Name	Recombinant Mouse F2RL1/PAR-2 Protein, N-GST & C-His
Host Species	E. coli
Purification	>90% as determined by SDS-PAGE.
Applications	ELISA, Immunogen, SDS-PAGE, WB, Bioactivity testing in progress
Other Names	GPR11, Thrombin receptor-like 1, G-protein coupled receptor 11, F2RL1, PAR-2, Coagulation factor II receptor-like 1, Proteinase-activated receptor 2, PAR2,
Accession No.	P55086
Target Species	Mus musculus (Mouse)
Calculated MW	40.42 kDa
Target Length	[Ser39-Gly73] & [Met214-Leu237] & [Ser350-Tyr399]
Formulation	Lyophilized from a solution in PBS pH 7.4, 0.02% NLS, 1mM EDTA, 4% Trehalose, 1% Mannitol.
Storage	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for one week. Store at -20 to -80°C for twelve months from the date of receipt.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.