

Recombinant Yellowfever mosquito Aed a 2/AeD7L1 Protein, C-His

Catalog No: #AP98233

Package Size: #AP98233-1 50ug #AP98233-2 100ug #AP98233-3 1mg

Description

Product Name	Recombinant Yellowfever mosquito Aed a 2/AeD7L1 Protein, C-His
Host Species	E. coli
Purification	>90% as determined by SDS-PAGE.
Applications	ELISA, Immunogen, SDS-PAGE, WB, Bioactivity testing in progress
Other Names	Long form salivary protein D7L1, AeD7, AeD7L1, 37 kDa salivary gland allergen Aed a 2, Protein D7, Aed a 2, D7, D7L1, AAEL006424, Aed a 2
Accession No.	P18153
Target Species	Aedes aegypti (Yellowfever mosquito) (Culex aegypti)
Calculated MW	36.19 kDa
Target Length	Ser18-Gln321
Tag Info	His-tag
Formulation	Lyophilized from a solution in PBS pH 7.4, 0.02% NLS, 1mM EDTA, 4% Trehalose, 1% Mannitol.
Storage	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for one week. Store at -20 to -80°C for twelve months from the date of receipt.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.