

Recombinant *Aspergillus fumigatus* Asp f 22/enoA Protein, N-His

Catalog No: #AP98236

Package Size: #AP98236-1 50ug #AP98236-2 100ug #AP98236-3 1mg

Description

Product Name	Recombinant <i>Aspergillus fumigatus</i> Asp f 22/enoA Protein, N-His
Host Species	Insect Cells
Purification	>90% as determined by SDS-PAGE.
Applications	ELISA, Immunogen, SDS-PAGE, WB, Bioactivity testing in progress
Other Names	Enolase, 4.2.1.11, 2-phospho-D-glycerate hydro-lyase, 2-phosphoglycerate dehydratase, Asp f 22, enoA, AFUA_6G06770
Accession No.	Q96X30
Target Species	<i>Aspergillus fumigatus</i> (strain ATCC MYA-4609 / CBS 101355 / FGSC A1100 / Af293) (<i>Neosartorya fumigata</i>)
Calculated MW	49.14 kDa
Target Length	Met1-Leu438(Val25Ala,Ala425Thr)
Tag Info	His-tag
Formulation	Lyophilized from a solution in PBS pH 7.4, 1mM EDTA, 4% Trehalose, 1% Mannitol.
Storage	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for one week. Store at -20 to -80°C for twelve months from the date of receipt.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.