

Recombinant American house dust mite Der f 1 Protein, C-His

Catalog No: #AP98238

Package Size: #AP98238-1 50ug #AP98238-2 100ug #AP98238-3 1mg

Description

Product Name	Recombinant American house dust mite Der f 1 Protein, C-His
Host Species	Yeast
Purification	>90% as determined by SDS-PAGE.
Applications	ELISA, Immunogen, SDS-PAGE, WB, Bioactivity testing in progress
Other Names	Der f 1, Der f 1, Der f 1 allergen
Accession No.	WWR30130.1
Target Species	Dermatophagoides farinae (American house dust mite)
Calculated MW	45.05 kDa
Target Length	Arg19-Met321
Tag Info	His-tag
Formulation	Lyophilized from a solution in PBS pH 7.4, 0.02% NLS, 1mM EDTA, 4% Trehalose, 1% Mannitol.
Storage	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for one week. Store at -20 to -80°C for twelve months from the date of receipt.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.