

Recombinant Peanut Ara h 1 Protein, N-His & C-Avi

Catalog No: #AP98240

Package Size: #AP98240-1 50ug #AP98240-2 100ug #AP98240-3 1mg

Description

Product Name	Recombinant Peanut Ara h 1 Protein, N-His & C-Avi
Host Species	E. coli
Purification	>90% as determined by SDS-PAGE.
Applications	ELISA, Immunogen, SDS-PAGE, WB, Bioactivity testing in progress
Other Names	Allergen Ara h 1, clone P41B, Allergen Ara h I, Ara h 1
Accession No.	P43238
Target Species	Arachis hypogaea (Peanut)
Calculated MW	49.69 kDa
Target Length	Ser170-Arg586
Formulation	Lyophilized from a solution in PBS pH 7.4, 0.02% NLS, 1mM EDTA, 4% Trehalose, 1% Mannitol.
Storage	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for one week. Store at -20 to -80°C for twelve months from the date of receipt.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.