

## Recombinant Cherry PRUA1/Pru av 1 Protein, C-His

Catalog No: #AP98249

Package Size: #AP98249-1 50ug #AP98249-2 100ug #AP98249-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

|                |                                                                                                                                                                       |
|----------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Product Name   | Recombinant Cherry PRUA1/Pru av 1 Protein, C-His                                                                                                                      |
| Host Species   | E. coli                                                                                                                                                               |
| Purification   | >90% as determined by SDS-PAGE.                                                                                                                                       |
| Applications   | ELISA, Immunogen, SDS-PAGE, WB, Bioactivity testing in progress                                                                                                       |
| Other Names    | Major allergen Pru av 1, Allergen Pru a 1, Pru av 1, PRUA1                                                                                                            |
| Accession No.  | O24248                                                                                                                                                                |
| Target Species | Prunus avium (Cherry) (Cerasus avium)                                                                                                                                 |
| Calculated MW  | 19.94 kDa                                                                                                                                                             |
| Target Length  | Met1-Asn160                                                                                                                                                           |
| Tag Info       | His-tag                                                                                                                                                               |
| Formulation    | Lyophilized from a solution in PBS pH 7.4, 0.02% NLS, 1mM EDTA, 4% Trehalose, 1% Mannitol.                                                                            |
| Storage        | Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for one week. Store at -20 to -80°C for twelve months from the date of receipt. |

Note: This product is for in vitro research use only and is not intended for use in humans or animals.