

Recombinant Honeybee Api m 4/Melittin Protein, N-His-SUMO

Catalog No: #AP98254

Package Size: #AP98254-1 50ug #AP98254-2 100ug #AP98254-3 1mg

Description

Product Name	Recombinant Honeybee Api m 4/Melittin Protein, N-His-SUMO
Host Species	E. coli
Purification	>90% as determined by SDS-PAGE.
Applications	ELISA, Immunogen, SDS-PAGE, WB, Bioactivity testing in progress
Other Names	Melittin, MEL, MLT, Allergen Api m 3, Allergen Api m III, Api m 4, MELT
Accession No.	P01501
Target Species	Apis mellifera (Honeybee)
Calculated MW	15.37 kDa
Target Length	Gly44-Gln69
Tag Info	His-tag
Formulation	Lyophilized from a solution in PBS pH 7.4, 0.02% NLS, 1mM EDTA, 4% Trehalose, 1% Mannitol.
Storage	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for one week. Store at -20 to -80°C for twelve months from the date of receipt.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.