

Recombinant Giant tiger prawn Pen m 1/TM1 Protein, C-His

Catalog No: #AP98255

Package Size: #AP98255-1 50ug #AP98255-2 100ug #AP98255-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Giant tiger prawn Pen m 1/TM1 Protein, C-His
Host Species	E. coli
Purification	>90% as determined by SDS-PAGE.
Applications	ELISA, Immunogen, SDS-PAGE, WB, Bioactivity testing in progress
Other Names	Tropomyosin, Major allergen Pen m 1, Tropomyosin, fast isoform, Tm-Penm-fast, Pen m 1.0101, TM1
Accession No.	A1KYZ2
Target Species	Penaeus monodon (Giant tiger prawn)
Calculated MW	33.90 kDa
Target Length	Met1-Tyr284
Tag Info	His-tag
Formulation	Lyophilized from a solution in PBS pH 7.4, 0.02% NLS, 1mM EDTA, 4% Trehalose, 1% Mannitol.
Storage	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for one week. Store at -20 to -80°C for twelve months from the date of receipt.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.