

Recombinant Mouse IgG2a Fc Fragment

Catalog No: #AP98266

Package Size: #AP98266-1 50ug #AP98266-2 100ug #AP98266-3 1mg

Description

Product Name	Recombinant Mouse IgG2a Fc Fragment
Host Species	Mammalian cells
Purification	>90% as determined by SDS-PAGE.
Applications	ELISA, Immunogen, SDS-PAGE, WB, Bioactivity testing in progress
Other Names	Ig gamma-2A chain C region, A allele, Immunoglobulin heavy chain gamma polypeptide, Ighg
Accession No.	P01863
Target Species	Mus musculus (Mouse)
Calculated MW	28.32 kDa
Target Length	Glu98-Lys330
Tag Info	Fc-tag
Formulation	Lyophilized from a solution in PBS pH 7.4, 1mM EDTA, 4% Trehalose, 1% Mannitol.
Storage	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for one week. Store at -20 to -80°C for twelve months from the date of receipt.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.