

Recombinant Human SCD Protein, N-His

Catalog No: #AP98267

Package Size: #AP98267-1 50ug #AP98267-2 100ug #AP98267-3 1mg

Description

Product Name	Recombinant Human SCD Protein, N-His
Host Species	E. coli
Purification	>90% as determined by SDS-PAGE.
Applications	ELISA, Immunogen, SDS-PAGE, WB, Bioactivity testing in progress
Other Names	Stearoyl-CoA desaturase, Fatty acid desaturase, Acyl-CoA desaturase, SCD, Delta(9)-desaturase, hSCD1, SCDOS, SCD1, FADS5, Delta-9 desaturase
Accession No.	O00767
Target Species	Homo sapiens (Human)
Calculated MW	33.26 kDa
Target Length	[Met1-Val72] & [Ala119-Tyr217] & [Val264-Gly359]
Tag Info	His-tag
Formulation	Lyophilized from a solution in PBS pH 7.4, 0.02% NLS, 1mM EDTA, 4% Trehalose, 1% Mannitol.
Storage	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for one week. Store at -20 to -80°C for twelve months from the date of receipt.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.