

Recombinant Human SFTPA1 Protein, N-His

Catalog No: #AP98269

Package Size: #AP98269-1 50ug #AP98269-2 100ug #AP98269-3 1mg

Description

Product Name	Recombinant Human SFTPA1 Protein, N-His
Host Species	E. coli
Purification	>90% as determined by SDS-PAGE.
Applications	ELISA, Immunogen, SDS-PAGE, WB, Bioactivity testing in progress
Other Names	Collectin-4, Pulmonary surfactant-associated protein A1, SFTP1, SP-A1, COLEC4, SP-A, PSP-A, SFTPA1, SFTPA1B, Alveolar proteinosis protein, SFTPA, PSPA, PSAP, 35 kDa pulmonary surfactant-associated protein
Accession No.	Q8IWL2
Target Species	Homo sapiens (Human)
Calculated MW	15.92 kDa
Target Length	Gln128-Phe248
Tag Info	His-tag
Formulation	Lyophilized from a solution in PBS pH 7.4, 0.02% NLS, 1mM EDTA, 4% Trehalose, 1% Mannitol.
Storage	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for one week. Store at -20 to -80°C for twelve months from the date of receipt.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.