

Biotin-labeled Recombinant Human IL6 Protein, C-His

Catalog No: #AP98270

Package Size: #AP98270-1 50ug #AP98270-2 100ug #AP98270-3 1mg

Description

Product Name	Biotin-labeled Recombinant Human IL6 Protein, C-His
Host Species	Mammalian cells
Purification	>90% as determined by SDS-PAGE.
Applications	ELISA, Immunogen, SDS-PAGE, WB, Bioactivity testing in progress
Other Names	Hybridoma growth factor, BSF-2, IFN-beta-2, CDF, IL6, Interleukin-6, IL-6, B-cell stimulatory factor 2, Interferon beta-2, IFNB2, CTL differentiation factor
Accession No.	P05231
Target Species	Homo sapiens (Human)
Calculated MW	25 kDa
Target Length	Val30-Met212
Tag Info	His-tag
Formulation	0.01M PBS, pH 7.4, 0.2% BSA, 50% Glycerol, 0.05% Proclin 300.
Storage	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for one week. Store at -20 to -80°C for twelve months from the date of receipt.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.