

# Recombinant Human RBFOX3 Protein, N-GST & C-His

Catalog No: #AP98275

Package Size: #AP98275-1 50ug #AP98275-2 100ug #AP98275-3 1mg

## Description

|                |                                                                                                                                                                       |
|----------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Product Name   | Recombinant Human RBFOX3 Protein, N-GST & C-His                                                                                                                       |
| Host Species   | E. coli                                                                                                                                                               |
| Purification   | >90% as determined by SDS-PAGE.                                                                                                                                       |
| Applications   | ELISA, Immunogen, SDS-PAGE, WB, Bioactivity testing in progress                                                                                                       |
| Other Names    | Neuronal nuclei antigen, RBFOX3, RNA binding protein fox-1 homolog 3, NeuN antigen, Fox-1 homolog C                                                                   |
| Accession No.  | A6NFN3                                                                                                                                                                |
| Target Species | Homo sapiens (Human)                                                                                                                                                  |
| Calculated MW  | 37.00 kDa                                                                                                                                                             |
| Target Length  | Pro98-Thr174                                                                                                                                                          |
| Formulation    | Lyophilized from a solution in PBS pH 7.4, 0.02% NLS, 1mM EDTA, 4% Trehalose, 1% Mannitol.                                                                            |
| Storage        | Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for one week. Store at -20 to -80°C for twelve months from the date of receipt. |

Note: This product is for in vitro research use only and is not intended for use in humans or animals.