

Recombinant Staphylococcus aureus entC2/SEC2 Protein, C-His

Catalog No: #AP98281

Package Size: #AP98281-1 50ug #AP98281-2 100ug #AP98281-3 1mg

Description

Product Name	Recombinant Staphylococcus aureus entC2/SEC2 Protein, C-His
Host Species	E. coli
Purification	>90% as determined by SDS-PAGE.
Applications	ELISA, Immunogen, SDS-PAGE, WB, Bioactivity testing in progress
Other Names	Enterotoxin type C-2, SEC2, entC2
Accession No.	P34071
Target Species	Staphylococcus aureus
Calculated MW	31.09 kDa
Target Length	Glu28-Gly266
Tag Info	His-tag
Formulation	Lyophilized from a solution in PBS pH 7.4, 0.02% NLS, 1mM EDTA, 4% Trehalose, 1% Mannitol.
Storage	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for one week. Store at -20 to -80°C for twelve months from the date of receipt.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.