

Recombinant HPV B19 VP1/Coat protein VP1 Protein, N-His

Catalog No: #AP98288

Package Size: #AP98288-1 50ug #AP98288-2 100ug #AP98288-3 1mg

Description

Product Name	Recombinant HPV B19 VP1/Coat protein VP1 Protein, N-His
Host Species	E. coli
Purification	>90% as determined by SDS-PAGE.
Applications	ELISA, Immunogen, SDS-PAGE, WB, Bioactivity testing in progress
Other Names	capsid protein 1, VP1
Accession No.	YP_004928146.1
Target Species	Human parvovirus B19 (HPV B19)
Calculated MW	18.26 kDa
Target Length	His80-Ser227
Tag Info	His-tag
Formulation	Lyophilized from a solution in PBS pH 7.4, 0.02% NLS, 1mM EDTA, 4% Trehalose, 1% Mannitol.
Storage	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for one week. Store at -20 to -80°C for twelve months from the date of receipt.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.