

Recombinant CA16 VP4/P1A Protein, N-GST & C-His

Catalog No: #AP98293

Package Size: #AP98293-1 50ug #AP98293-2 100ug #AP98293-3 1mg

Description

Product Name	Recombinant CA16 VP4/P1A Protein, N-GST & C-His
Host Species	E. coli
Purification	>90% as determined by SDS-PAGE.
Applications	ELISA, Immunogen, SDS-PAGE, WB, Bioactivity testing in progress
Other Names	Genome polyprotein, Capsid protein VP4, P1A, Virion protein 4
Accession No.	Q9QF31
Target Species	Coxsackievirus A16 (strain Tainan/5079/98)
Calculated MW	35.59 kDa
Target Length	Gly2-Lys69
Formulation	Lyophilized from a solution in PBS pH 7.4, 0.02% NLS, 1mM EDTA, 4% Trehalose, 1% Mannitol.
Storage	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for one week. Store at -20 to -80°C for twelve months from the date of receipt.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.