

Recombinant CV-A10 VP1/P1D Protein, N-GST/C-His

Catalog No: #AP98296

Package Size: #AP98296-1 50ug #AP98296-2 100ug #AP98296-3 1mg

Description

Product Name	Recombinant CV-A10 VP1/P1D Protein, N-GST/C-His
Host Species	E. coli
Purification	>90% as determined by SDS-PAGE.
Applications	ELISA, Immunogen, SDS-PAGE, WB, Bioactivity testing in progress
Other Names	Genome polyprotein, Capsid protein VP1, P1D, Virion protein 1
Accession No.	USC32119.1
Target Species	Coxsackievirus A10 (CV-A10)
Calculated MW	49.07 kDa
Target Length	Thr78-Val261
Tag Info	GST-tag
Formulation	Lyophilized from a solution in PBS pH 7.4, 0.02% NLS, 1mM EDTA, 4% Trehalose, 1% Mannitol.
Storage	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for one week. Store at -20 to -80°C for twelve months from the date of receipt.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.