

Recombinant EV68 P3C/Protease 3C Protein, N-His

Catalog No: #AP98297

Package Size: #AP98297-1 50ug #AP98297-2 100ug #AP98297-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant EV68 P3C/Protease 3C Protein, N-His
Host Species	E. coli
Purification	>90% as determined by SDS-PAGE.
Applications	ELISA, Immunogen, SDS-PAGE, WB, Bioactivity testing in progress
Other Names	Genome polyprotein, Protease 3C, 3.4.22.28, Picornain 3C, P3C
Accession No.	Q68T42
Target Species	Human enterovirus D68 (EV68) (EV-68)
Calculated MW	22.50 kDa
Target Length	Gly1549-Gln1731
Tag Info	His-tag
Formulation	Lyophilized from a solution in PBS pH 7.4, 0.02% NLS, 1mM EDTA, 4% Trehalose, 1% Mannitol.
Storage	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for one week. Store at -20 to -80°C for twelve months from the date of receipt.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.