

Recombinant EV68 VP2/P1B Protein, N-His

Catalog No: #AP98300

Package Size: #AP98300-1 50ug #AP98300-2 100ug #AP98300-3 1mg

Description

Product Name	Recombinant EV68 VP2/P1B Protein, N-His
Host Species	E. coli
Purification	>90% as determined by SDS-PAGE.
Applications	ELISA, Immunogen, SDS-PAGE, WB, Bioactivity testing in progress
Other Names	Genome polyprotein, Capsid protein VP2, P1B, Virion protein 2
Accession No.	Q68T42
Target Species	Human enterovirus D68 (EV68) (EV-68)
Calculated MW	29.97 kDa
Target Length	Ser70-Gln317
Tag Info	His-tag
Formulation	Lyophilized from a solution in PBS pH 7.4, 0.02% NLS, 1mM EDTA, 4% Trehalose, 1% Mannitol.
Storage	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for one week. Store at -20 to -80°C for twelve months from the date of receipt.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.