

Recombinant EV71 P2A/Protease 2A Protein, N-His

Catalog No: #AP98301

Package Size: #AP98301-1 50ug #AP98301-2 100ug #AP98301-3 1mg

Description

Product Name	Recombinant EV71 P2A/Protease 2A Protein, N-His
Host Species	E. coli
Purification	>90% as determined by SDS-PAGE.
Applications	ELISA, Immunogen, SDS-PAGE, WB, Bioactivity testing in progress
Other Names	Genome polyprotein, Protease 2A, P2A, 3.4.22.29, Picornain 2A, Protein 2A
Accession No.	Q66479
Target Species	Human enterovirus 71 (strain 7423/MS/87) (EV71) (EV-71)
Calculated MW	18.75 kDa
Target Length	Gly863-Gln1012
Tag Info	His-tag
Formulation	Lyophilized from a solution in PBS pH 7.4, 0.02% NLS, 1mM EDTA, 4% Trehalose, 1% Mannitol.
Storage	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for one week. Store at -20 to -80°C for twelve months from the date of receipt.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.