

Recombinant EV71 VP1/P1D Protein, N-Fc

Catalog No: #AP98305

Package Size: #AP98305-1 50ug #AP98305-2 100ug #AP98305-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant EV71 VP1/P1D Protein, N-Fc
Host Species	E. coli
Purification	>90% as determined by SDS-PAGE.
Applications	ELISA, Immunogen, SDS-PAGE, WB, Bioactivity testing in progress
Other Names	Genome polyprotein, Capsid protein VP1, P1D, Virion protein 1
Accession No.	F5AU45
Target Species	Enterovirus A (GZ2010-156)
Calculated MW	60.26 kDa
Target Length	Gly566-Leu862
Tag Info	Fc-tag
Formulation	Lyophilized from a solution in PBS pH 7.4, 0.02% NLS, 1mM EDTA, 4% Trehalose, 1% Mannitol.
Storage	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for one week. Store at -20 to -80°C for twelve months from the date of receipt.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.