

Recombinant KSI tag Control Protein (N-His)

Catalog No: #AP98308

Package Size: #AP98308-1 50ug #AP98308-2 100ug #AP98308-3 1mg

Description

Product Name	Recombinant KSI tag Control Protein (N-His)
Host Species	E. coli
Purification	>90% as determined by SDS-PAGE.
Applications	ELISA, Immunogen, SDS-PAGE, WB, Bioactivity testing in progress
Other Names	KSI tag, KSI-tag, Steroid Delta-isomerase, Delta(5)-3-ketosteroid isomerase
Accession No.	P00947
Target Species	Comamonas testosteroni (Pseudomonas testosteroni)
Calculated MW	18.04 kDa
Target Length	Met1-Ala123
Tag Info	His-tag
Formulation	Lyophilized from a solution in PBS pH 7.4, 0.02% NLS, 1mM EDTA, 4% Trehalose, 1% Mannitol.
Storage	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for one week. Store at -20 to -80°C for twelve months from the date of receipt.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.