

Recombinant Influenza A virus (H10N5) NA/Neuraminidase Protein, C-His

Catalog No: #AP98322

Package Size: #AP98322-1 50ug #AP98322-2 100ug #AP98322-3 1mg

Description

Product Name	Recombinant Influenza A virus (H10N5) NA/Neuraminidase Protein, C-His
Host Species	Mammalian cells
Purification	>90% as determined by SDS-PAGE.
Applications	ELISA, Immunogen, SDS-PAGE, WB, Bioactivity testing in progress
Other Names	neuraminidase, NA
Accession No.	AYD76504.1
Target Species	Influenza A virus (A/wild birds/Hubei/102/2014 (H10N5))
Calculated MW	51.28 kDa
Target Length	Ser35-Met472
Tag Info	His-tag
Formulation	Lyophilized from a solution in PBS pH 7.4, 1mM EDTA, 4% Trehalose, 1% Mannitol.
Storage	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for one week. Store at -20 to -80°C for twelve months from the date of receipt.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.