

# Recombinant Influenza A virus (H10N7) HA/Hemagglutinin Protein, C-His

Catalog No: #AP98325

Package Size: #AP98325-1 50ug #AP98325-2 100ug #AP98325-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

|                |                                                                                                                                                                       |
|----------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Product Name   | Recombinant Influenza A virus (H10N7) HA/Hemagglutinin Protein, C-His                                                                                                 |
| Host Species   | Mammalian cells                                                                                                                                                       |
| Purification   | >90% as determined by SDS-PAGE.                                                                                                                                       |
| Applications   | ELISA, Immunogen, SDS-PAGE, WB, Bioactivity testing in progress                                                                                                       |
| Other Names    | Hemagglutinin, Hemagglutinin HA1 chain, Hemagglutinin HA2 chain, HA                                                                                                   |
| Accession No.  | Q0A448                                                                                                                                                                |
| Target Species | Influenza A virus (strain A/Duck/Germany/1949 (H10N7))                                                                                                                |
| Calculated MW  | 63.78 kDa                                                                                                                                                             |
| Target Length  | Leu17-Asp524                                                                                                                                                          |
| Tag Info       | His-tag                                                                                                                                                               |
| Formulation    | Lyophilized from a solution in PBS pH 7.4, 1mM EDTA, 4% Trehalose, 1% Mannitol.                                                                                       |
| Storage        | Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for one week. Store at -20 to -80°C for twelve months from the date of receipt. |

Note: This product is for in vitro research use only and is not intended for use in humans or animals.