

Recombinant Influenza A virus (H14N6) HA/Hemagglutinin Protein, C-His

Catalog No: #AP98334

Package Size: #AP98334-1 50ug #AP98334-2 100ug #AP98334-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Influenza A virus (H14N6) HA/Hemagglutinin Protein, C-His
Host Species	Mammalian cells
Purification	>90% as determined by SDS-PAGE.
Applications	ELISA, Immunogen, SDS-PAGE, WB, Bioactivity testing in progress
Other Names	hemagglutinin, HA
Accession No.	AGE03043.1
Target Species	Influenza A virus (A/mallard/Wisconsin/10OS3941/2010 (H14N6))
Calculated MW	64.73 kDa
Target Length	Ser17-Asp531
Tag Info	His-tag
Formulation	Lyophilized from a solution in PBS pH 7.4, 1mM EDTA, 4% Trehalose, 1% Mannitol.
Storage	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for one week. Store at -20 to -80°C for twelve months from the date of receipt.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.