

Recombinant Influenza A virus (H17N10) HA/Hemagglutinin Protein, N-His

Catalog No: #AP98337

Package Size: #AP98337-1 50ug #AP98337-2 100ug #AP98337-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Influenza A virus (H17N10) HA/Hemagglutinin Protein, N-His
Host Species	E. coli
Purification	>90% as determined by SDS-PAGE.
Applications	ELISA, Immunogen, SDS-PAGE, WB, Bioactivity testing in progress
Other Names	hemagglutinin, HA
Accession No.	AFC35418.1
Target Species	Influenza A virus (A/little yellow-shouldered bat/Guatemala/153/2009 (H17N10))
Calculated MW	29.00 kDa
Target Length	Gly57-Lys293
Tag Info	His-tag
Formulation	Lyophilized from a solution in PBS pH 7.4, 0.02% NLS, 1mM EDTA, 4% Trehalose, 1% Mannitol.
Storage	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for one week. Store at -20 to -80°C for twelve months from the date of receipt.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.